



## Data Collection Worksheet

**Please Note:** The Data Collection Worksheet (DCW) is a tool to aid integration of a PhenX protocol into a study. The PhenX DCW is not designed to be a data collection instrument. Investigators will need to decide the best way to collect data for the PhenX protocol in their study. Variables captured in the DCW, along with variable names and unique PhenX variable identifiers, are included in the PhenX Data Dictionary (DD) files.

TA: 6:03

PAT: 2

Voxel size: 1.0 x 1.0 x 1.0 mm

Rel. SNR: 1.00

### Properties

Prio Recon .....Off  
Before  
measurement  
After  
measurement  
Load to viewer .....On  
Inline movie .....Off  
Auto store images .....On  
Load to stamp .....Off  
segments  
Load images to .....Off  
graphic segments  
Auto open inline .....Off  
display  
Load to stamp .....Off  
segments  
Load images to .....Off  
graphic segments  
Auto open inline .....Off  
display  
Start .....On  
measurement  
without further

preparation  
 Wait for user to .....On  
 start  
 Start .....single  
 measurements

## Routine

Slab group 1  
 Slab .....1  
 Dist. ....50%  
 Factor .....L0.0  
 Position A14.2 H0.7  
 Orientation .....S > C0.2  
 Phase enc. ....A >> P  
 dir. ....0.00 deg  
 Rotation  
 Phase .....0%  
 oversampling  
 Slice .....0.0%  
 oversampling  
 Slices per slab .....176  
 FoV read .....256 mm  
 FoV phase .....100.0%  
 Slice thickness .....1.00 mm  
 TR .....2530 ms  
 TE 1 .....3.31 ms  
 TE 2 .....6.99 ms  
 Averages .....1  
 Concatenations .....1  
 Filter .....Prescan  
 Normalize  
 Coil elements .....HEA;HEP

## Contrast

Magn. ....Non-sel.  
 preparation IR  
 TI .....1100 ms  
 Flip angle .....7.0 deg  
 Fat suppr. ....None

Water suppr. ....None  
 Averaging .....Long term  
 mode  
 Reconstruction .....Magnitude  
 Measurements .....1  
 Multiple series .....Each  
 measurement

## Resolution

Base resolution .....256  
 Phase .....100%  
 resolution  
 Slice resolution .....100%  
 Phase partial .....Off  
 Fourier  
 Slice partial ..... Off  
 Fourier  
 Interpolation .....Off  
 PAT mode .....GRAPPA  
 Accel. ....2  
 factor PE  
 Ref. lines .....32  
 PE  
 Matrix Coil .....Auto  
 Mode (Triple)  
 Reference .....Integrated  
 scan mode  
 Image Filter .....Off  
 Distortion .....Off  
 Corr.  
 Unfiltered .....Off  
 images  
 Prescan .....On  
 Normalize  
 Normalize .....Off  
 B1 filter .....Off  
 Raw filter .....Off  
 Elliptical filter .....Off

## Geometry

Multi-slice .....Single shot  
mode  
Series .....Interleaved

**System**

Body .....Off  
HEP .....On  
HEA .....On  
Positioning .....FIX  
mode  
Table .....H  
position  
Table .....0 mm  
position  
MSMA .....S - C - T  
Sagittal .....R >> L  
Coronal .....A >> P  
Transversal .....F >> H  
Save .....Off  
uncombined  
Coil .....Adaptive  
Combine Combine  
Mode  
AutoAlign .....Head >  
Brain  
Auto Coil .....Default  
Select  
Shim mode .....Standard  
Adjust with .....Off  
body coil  
Confirm .....Off  
freq.  
adjustment  
Assume .....Off  
Silicone  
? Ref. .....0.000 V  
amplitude 1H  
Adjustment .....Auto  
Tolerance  
Adjust  
volume  
Position .....L0.0

A14.2 H0.7

Orientation .....S > C0.2  
Rotation .....0.00 deg  
F >> H .....256 mm  
A >> P .....256 mm  
R >> L .....176 mm

Physio

1st .....None  
Signal/Mode  
Dark blood .....Off

Inline

Subtract .....Off  
Std-Dev-Sag .....Off  
Std-Dev-Cor .....Off  
Std-Dev-Tra .....Off  
Std-Dev-Time .....Off  
MIP-Sag .....Off  
MIP-Cor .....Off  
MIP-Tra .....Off  
MIP-Time .....Off  
Save original  
images .....On

Sequence

Introduction .....On  
Dimension .....3D  
Elliptical  
scanning .....Off  
Asymmetric  
echo .....Off  
Contrasts .....2  
Bandwidth 1 .....195  
Hz/Px  
Bandwidth 2 .....651

	Hz/Px
Flow comp. 1	No
Flow comp. 2	No
Echo spacing	10.1 ms
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Add. scale factor	8.0
Gradient spoiling	Integral
Gradient moment factor	3.0
Siemens reconstruction	On
Save raw k-space data	Off
Averaging	RMS

Protocol source: <https://www.phenxtoolkit.org/protocols/view/660501>